






























XPZ LABORATORINIŲ INDŲ PLOVIMO MAŠINOS IR PRIEDAI





Eil. Nr.	Modelis	Pavadinimas	Vaizdas	Aprašymas	Techninės charakteristikos
1.	Aurora-1	Glassware washer		Use for cleaning and dry the glass, ceramics, metal, plastic laboratory glassware, Also for Flask, Vials, wide-necked bottle, triangle flask, measuring flask, colorimetric tube, culture dish ect.	1. Working environment: 5-40°C 2. Adaptor: 220V/380V, Single Phase voltage 50Hz 3. Circulating pump: 900W 4. Circulating system: 0-400L/min 5. Temperature sensor: PT1000 dual detector 6. Heating power: 2000W/4000W 7. Housing&Liner: 316L stainless steel 8. Size: 880*600*640mm 9. LCD color: Black/Orange/Blue/Green 10. Packing: Carton box
2.	Aurora-F1	Glassware washer		Use for cleaning and dry the glass, ceramics, metal, plastic laboratory glassware, Also for Flask, Vials, wide-necked bottle, triangle flask, measuring flask, colorimetric tube, culture dish ect.	1. Working environment: 5-40°C 2. Adaptor: 220V/380V, Single Phase voltage 50Hz 3. Circulating pump: 900W 4. Circulating system: 0-400L/min 5. Temperature sensor: PT1000 dual detector 6. Heating power: 2000W/4000W 7. Housing&Liner: 316L stainless steel 8. Hot dry system: HEPA 9. Size: 880*920*650mm 10. LCD color: Black/Orange/Blue/Green 11. Packing: Carton box
3.	DZ-101	Full injection cleaning rack		Use for flask, vials, test tube, triangle flask, pipette, round-bottle flask ect.	Max cleaning quantity: vials/test tube \geq 72 , with fix claw \geq 28.
4.	DZ-102	Full injection cleaning rack		Use for flask, vials, test tube, triangle flask, pipette, round-bottle flask ect.	Max cleaning quantity: vials/test tube \geq 72 , with fix claw \geq 28. Free solid wash powder nozzle
5.	DZ-103	Full injection cleaning rack		Use for flask, colorimetric tube, test tube ect.	Max cleaning quantity: 40-52 , Free solid wash powder nozzle

6.	DZ-104	Full injection cleaning rack		Use for flask,colorimetric tube,test tube ect.	Max cleaning quantity: 40-52 , Free solid wash powder nozzle
7.	DZ-105	Full injection cleaning rack		Use for flask,colorimetric tube,triangle flask,blue cap bottle,separating funnel,test tube ect.	Max cleaning quantity: 34-50 , Free solid wash powder nozzle
8.	DZ-106	Full injection cleaning rack		Use for flask,colorimetric tube,triangle flask,blue cap bottle,separating funnel,test tube ect.	Max cleaning quantity: 34-50 , Free solid wash powder nozzle
9.	DZ-107	Semi injection cleaning rack		Use for flask,colorimetric tube,triangle flask,blue cap bottle,separating funnel,test tube beaker,measuring flask,wide-necked bottle ect	Max cleaning quantity: 17-25
10.	DZ-108	Semi injection cleaning rack		Use for flask,colorimetric tube,triangle flask,blue cap bottle,separating funnel,test tube beaker,measuring flask,wide-necked bottle ect	Max cleaning quantity: 17-25
11.	DZ-109	Semi injection cleaning rack		Use for flask,colorimetric tube,, separating funnel,test tube beaker,measuring flask,wide-necked bottle ect	Max cleaning quantity: 20-26 , Free solid wash powder nozzle
12.	DZ-110	Semi injection cleaning rack		Use for flask,colorimetric tube,, separating funnel,test tube beaker,measuring flask,wide-necked bottle ect	Max cleaning quantity: 20-26 , Free solid wash powder nozzle
13.	DZ-111	Full injection cleaning rack		Use for pipette	Max cleaning quantity: 116 Pipette max height:within 460mm

14.	DZ-112	Full injection cleaning rack		Use for pipette	Max cleaning quantity: 116 max height:within 460mm	Pipette
15.	DZ-113	Full injection cleaning rack		Use for absorption bottle which is for environmental monitoring	Max cleaning quantity: 20-40	
16.	DZ-114	Full injection cleaning rack		Use for absorption bottle which is for environmental monitoring	Max cleaning quantity: 20-40	
17.	SZ-101	Upper injection cleaning rack		Use for Flask(less than 250ml), triangle flask, round-bottle flask, colorimetric tube ect.	Max cleaning quantity: 34	
18.	SZ-102	Upper injection cleaning rack		Use for Flask(less than 250ml), triangle flask, round-bottle flask, colorimetric tube ect.	Max cleaning quantity: 34	
19.	SZ-103	Upper injection cleaning rack		Use for Vials, centrifuge tube ect.	Max cleaning quantity: 96-230	
20.	SZ-104	Upper injection cleaning rack		Use for Vials, centrifuge tube ect.	Max cleaning quantity: 96-230	
21.	B-101	Fixed Claw		Fix the glassware	Fix the glassware	

22.	G-101	Stainless steel net (for vials)		special design for preventing the vial move out.	special design for preventing the vial move out.
23.	G-201	Stainless steel net (half)		special design for preventing the glassware move out.	special design for preventing the glassware move out.
24.	G-301	Stainless steel net (one fourth)		special design for preventing the glassware move out	special design for preventing the glassware move out
25.	ST-101	tray(bottom)		sitable for different frame	sitable for different frame
26.	TSXF-101	frame(Semi-rimless)		For beaker, wide-necked bottle, funnel, glass rod, cap (with handle)	For beaker, wide-necked bottle, funnel, glass rod, cap (with handle)
27.	XZ-101	Full injection cleaning rack(bottom)		Use for Flask (less than 250ml), triangle flask, round-bottle flask, colorimetric tube ect.	Max cleaning quantity: 34
28.	XZ-103	Full injection cleaning rack(bottom)		Use for Vials, centrifuge tube ect.	Max cleaning quantity: 96-230
29.	TSXF-204	frame(Full-rim)		for small glassware	for small glassware

30.	TXF-105	frame(Semi-rimless)		for small glassware	for small glassware
31.	TSXF-401	frame(one fourth)		For test tube, colorimetric tube ect.	Max cleaning quantity: ≥ 160
32.	TXF-101	frame(Full-rim)		For culture dish	Max cleaning quantity: ≥ 56
33.	TXF-102	frame(Full-rim)		For culture dish (only use with TXF-101)	Max cleaning quantity: ≥ 56
34.	TXF-103	frame(Botom, full-rim)		For reagent bottle, blue cap bottle ect.	Max cleaning quantity: ≥ 9
35.	TXF-104	frame(Botom, full-rim)		For reagent bottle, blue cap bottle ect.	Max cleaning quantity: ≥ 20
36.	Z-101	pipette accessories		Only use for pipette	Only use for pipette
37.	TXF-201	Semi cleaning rack		For beaker, measuring flask, triangle flask ect.	pinch cock:28

38.	TXF-202	Semi cleaning rack		For beaker, measuring flask, triangle flask ect.	pinch cock:28
39.	TXF-203	Semi cleaning rack		For beaker, measuring flask ect.	Max cleaning quantity:21
40.	TSXF-202	Semi cleaning rack		For watch glass	Max cleaning quantity: ≥ 105
41.	TSXF-203	Semi cleaning rack		For watch glass	Max cleaning quantity: ≥ 44