

Whitening products

Professional whitening



Tecno-Gaz products

The products used have been tested and certified over time according to the corresponding European standards. **Tecno-Gaz has chosen a strategy of controlling the prices of its consumer products, allowing all professionals to be able to offer the treatments to a greater number of clients.** Undoubtedly an important and strategic choice that will ensure that every dentist has the chance to make whitening an important business aspect for their surgery, ensuring client satisfaction and loyalty.



BasicSet

- Beta: n°5 blisters containing 38% liquid hydrogen peroxide
- Alfa: n° 5 syringes containing 0.75g of activator
- LiquidMask: n°3 syringes with 3.5g with liquid dam for polymerisation
- 1 Mixing spatula;
- 10 Micro brushes for application;
- 10 flow tips for dam

Art. 4418S2




Desens

Desensitiser kit potassium nitrate

Includes 10 single patient use bottles 0.5 ml each + 10 brushes.

Art. 4412S2



As ergonomic as possible

- Below the head there is a master switch that switches the machine on and off and starts and stops the treatment.
- Switching the machine on will activate a timed session, lasting 10 minutes, which will then automatically switch off.

As practical as possible

Everything is contained in a practical plastic case useful for transporting the device, from one surgery to another.

Clip

Weight	5 Kg
Single-phase power supply	230 Volt
Frequency	50 Hz
Maximum absorbed power	50 VA
Fuse	F 1 A

Light emission characteristics

Radiation wavelength	495 nm
Maximum radiated light power	30.000 mw
Medical device	Class I 93/42/CEE
Electric safety	Class I Type B CEI EN 60601-1
Type of protection	PE

