

# MediCanvas

## TFT-LCD MEDiCanvas X-ray view boxes

Tecno-Gaz MediCanvas X-ray view boxes use TFT-LCD technology and are at the top of the market as regards quality, aesthetics and technology.

A prestigious product that will not only improve your operating efficiency, giving your surgery added value, but that also has a striking aesthetic impact.

All MediCanvas X-ray view boxes are already fitted out for wallmounting, and items 7072-N 7074-N may be used standing on the table, with the aid of the appropriate support. These negatoscopes are equipped with a device that allows automatic FAS turning on and turning off when you insert the x Ray.

## Advantages



### Excellent light uniformity

Computer screen technology is used to provide total light **uniformity** and this facilitates the reading of the X-ray and avoids the operator straining their eyes.



### Easy X-ray reading

MediCanvas uses CCFL technology that **guarantees the emission of cold blue light at 8.600° Kelvin**, allowing for better definition when consulting X-rays.



### Reduced eye fatigue

The CCFL system uses a frequency of 40,000 times per second. **This method avoids any movement of the lamp**, unlike in traditional X-ray viewing boxes, **noticeably reducing eye fatigue**.



### Long-life lamp

The MediCanvas lighting system **lasts for approximately 10 years without specific maintenance**, thanks to the CCFL system that **ensures time and cost savings**.



### MediCanvas design and quality

**MediCanvas X-ray viewing boxes are available in 4 different sizes, designed for every need.** The main distinguishing feature of these products is their reduced thickness that varies from 24 to 39 mm. Their absolutely extraordinary aesthetic impact and their presence provides a very professional look to any environment.



### Ask the expert

#### Who needs this?

All dental surgeries.

#### Why MediCanvas X-ray viewing boxes?

An innovative product that makes use of technology to improve your work. A pleasantly designed, high-quality product.



Frame and click!


Use your mobile phone to get further information about this issue.



Art. 7072-N



Art. 7074-N



Art. 7076-N



Art. 7078-N



## Optional

### Frame 7077-N

**Counter support** for  
MediCanvas X-ray viewing box  
Art. 7072-N

### Frame 7079-N

**Counter support** for  
MediCanvas X-ray viewing box  
Art. 7074-N



## X-ray viewing box MediCanvas with TFT-LCD technology

### Frame 7072-N

Width	498 mm
Height	482 mm
Depth	28 mm

### Viewing Area 14x17"

Width	356 mm
Height	423 mm

### Frame 7074-N

Width	854 mm
Height	482 mm
Depth	28 mm

### Viewing Area 28x17"

Width	711 mm
Height	423 mm

### Frame 7076-N

Width	1230 mm
Height	482 mm
Depth	39 mm

### Viewing Area 43x17"

Width	1087 mm
Height	423 mm

### Frame 7078-N

Width	1606 mm
Height	482 mm
Depth	39 mm

### Viewing Area 58x17"

Width	1463 mm
Height	423 mm